# Hotgen

## Coronavirus (2019-nCoV)-Antigen test-Instructions for Use

**REF**  
HGGC134S0101 (17/KR)  
HGGC134S0105 (57/KR)  
HGGC134S0120 (207/KR)  
HGGC134S0140 (407/KR)

**SPECIMEN COLLECTION**

Children and young people under the age of 18 should be supported by an adult.

### COMPONENTS

Instructions for Use SARS-CoV-2 antigen test cassette


### TEST STEPS

**1** Please wash and dry your hands thoroughly before the test.

**8** Soak the swab in the buffer for at least 15 seconds, stir the swab several times, and squeeze the swab head 3 times.

**10** Close the sampling tube with the tube lid

**11** Open the foil pouch and place the test cassette on a flat surface.

**12** Open the tube lid at the front end of the sampling tube, and drop 4 drops of the treated sample into the sample well (S) of the test cassette.

**15 min** After reacting at room temperature for 15 minutes, observe the results. A result after 30 minutes is invalid.

**2** Read the instructions for use carefully.

**3** Please scan the QR code to watch the operation video.

**4** Take the swab out of the package and do not touch the sampling end.

**5** Carefully insert the swab 1.5 cm into the nostril until resistance is felt. Under moderate pressure, swab the surface of the nostril 4 to 6 times for at least 15 seconds.

**6** Repeat sampling with the same swab in the other nostril.

**7** Open the larger end of the sampling tube, and insert the swab after collecting the sample into the sampling tube.

**14** After the test, pack all the components of this test in the Biohazard specimen bag for contaminated waste and dispose of this bag closed with the residual waste. Not reusable.

**15** Wash or disinfect hands again.

**10** Close the sampling tube with the tube lid

**11** Open the foil pouch and place the test cassette on a flat surface.

**12** Open the tube lid at the front end of the sampling tube, and drop 4 drops of the treated sample into the sample well (S) of the test cassette.

**15 min** After reacting at room temperature for 15 minutes, observe the results. A result after 30 minutes is invalid.

### INTENDED USE

This kit for self-testing is used for in vitro qualitative determination of SARS-CoV-2 antigen in human anterior nasal swab samples. It can be used for rapid investigation of SARS-CoV-2 antigen from individuals suspected of COVID-19 within the first seven days of the onset of symptoms and the individuals without symptoms. Children and young people under the age of 18 should be supported by an adult.

**MATERIALS PROVIDED**

Components	17/KR	57/KR	207/KR	407/KR
SARS-CoV-2 antigen test cassette	1 test	5 tests	20 tests	40 tests
Sample extraction buffer	1 pcs	5 pcs	20 pcs	40 pcs
Disposable Sterile Swab	1 test	5 tests	20 tests	40 tests
Biohazard specimen bag	1 pcs	5 pcs	20 pcs	40 pcs
Instructions for use	1 pcs	1 pcs	1 pcs	1 pcs

### PERFORMANCE

**1. Exogenous/Endogenous Interference Substances studies:** There was no interference for potential interfering substances listed below. (1)Exogenous factor

No.	Exogenous factor	Interfering substances	Test conc.
1		Phenylephrine	100 µg/mL
2		Dexamethasone	10 µg/mL
3		Saline Nasal Spray 0.9%	10 µg/mL
4		Dexamethasone	2 µg/mL
5		Flunisolide	0.2 µg/mL
6		Tramancolone acetate	0.2 µg/mL
7		Mometasone	0.5 µg/mL

**(2)Endogenous factor**

No.	Exogenous factor	Interfering substances	Test conc.
1	Autimmune disease	Human anti-mouse antibody (HAMA)	800 ng/mL
2	Serum protein	Whole Blood (human), EDTA anticoagulated	10% (v/v)

**2. Cross-Reactivity & Microbial interference:** There was no cross-reaction and interference with the potential cross-reacting microorganisms listed below.

No.	Crossing reacting substance	Strain	Concentration of cross reacting substance
1		HKU1	2 × 10 <sup>7</sup> TCID50/mL
2		Z99E	2 × 10 <sup>7</sup> TCID50/mL
3		OC43	2 × 10 <sup>7</sup> TCID50/mL
4		NL63	2 × 10 <sup>7</sup> TCID50/mL
5		SARS	2 × 10 <sup>7</sup> TCID50/mL
6		MERS	2 × 10 <sup>7</sup> TCID50/mL
7		H1N1	2 × 10 <sup>7</sup> TCID50/mL
8		H3N2	2 × 10 <sup>7</sup> TCID50/mL
9		H5N1	2 × 10 <sup>7</sup> TCID50/mL
10		H7N9	2 × 10 <sup>7</sup> TCID50/mL

**EXPLANATION OF SYMBOL**

Symbol	Meaning
Use by date	Consult instruction for use
Content Sufficient For re-tests	Temperature limitation
Manufacturing date	Caution
CE Marking - IVD 98/79/EC	Authorized representative in the European Community
For In Vitro Diagnostic Use	Keep away from sunlight
Do not use if package is damaged	Keep dry
CE mark	Sterilized using ethylene oxide

### ACCESSORIES

Accessory	Manufacturer	EC-Representative	CE Mark
Swab	Zhejiang Gongdong Medical Technology Co., Ltd. NO.10 Binyuan Ave.Huang yan 30020 Taizhou, Zhejiang PEOPLES REPUBLIC OF CHINA	Shanghai International Holding corp GmbH (Europe) Eiffelstraße 80 20537 Hamburg Germany	CE

**LITERATURE REFERENCE**

1. Wu C, Liu Y, Yang Y, Zhang P, Zhong W, Wang Y, et al. (February 2020). "Analysis of therapeutic targets for SARS-CoV-2 and discovery of potential drugs by computational methods". Acta Pharmaceutica Sinica B. doi:10.1016.

### FREQUENTLY ASKED QUESTIONS (FAQ)

- When can I test myself?** You can always test yourself whether you have symptoms or not. Please note that the test result is a snapshot that is valid for this point in time. Tests should therefore be repeated in accordance with the requirements of the responsible authorities.
- What do I have to pay attention to in order to get the most exact test result possible?** Always follow the instructions for use exactly. Perform the test immediately after collecting the sample. Put the drops from the sample tube only into the designated well of the test cassette. Dispense four drops from the sample tube. Too many or too few drops can lead to an incorrect or invalid test result.
- How does the test work?** The H-protein of the SARS-CoV-2 virus reacts with the coating of the test line and leads to a color change, i.e. a red line appears. If the sample does not contain any virus proteins or antigens, no red test line (T) appears.
- The test strip is very discolored. Why is this or what am I doing wrong?** The reason for that too large a quantity of drops has been dispensed from the test tube into the test cassette well. The indicator strip can only hold a limited amount of liquid. If the control line does not appear or the test strip is severely discolored, please repeat the test with a new test kit according to the instructions for use.
- My result is positive. What should I do?** If a horizontal colored line is visible in the control area (C) as well as in the test area (T), your result is positive. There is currently a suspicion of a COVID-19 infection. Immediately take self-isolation measures in accordance with the local guidelines and immediately contact your general practitioner / doctor or the local health department in accordance with the instructions of your local authorities. Your test result will be checked by a PCR confirmation test and the next steps will be explained to you.
- My result is negative. What should I do?** If only a horizontal colored line is visible in the control area (C), this can mean that you are negative or that the viral load is too low to be recognized by the test. If you experience symptoms such as headaches, migraines, fever, loss of sense of smell and taste, contact the nearest medical facility using the rules of your local authority. In addition, you can repeat the test with a new test kit. In case of suspicion, repeat the test after 1-2 days, as the coronavirus cannot be precisely detected in all phases of an infection. Distance and hygiene rules must still be observed.
- Can this test cassette be reused or used by more than one person?** This test cassette is for single use and cannot be reused or used by more than one person.
- Where do I dispose of the product?** The test kit can be disposed of with normal household waste.

3. Coronavirus (2019-nCoV)-Antigen test - was compared to the RT-PCR test results. The clinical study involved 400 samples.

Item	Positive samples	Negative samples
Number	131	269
Sensitivity	96.95%	
Specificity	98.88%	
Accuracy	98.25%	

The layperson study involved 203 samples in total, including 24 positive samples and 179 negative samples. The feasibility study demonstrated that: 98.03% of non-professionals carried out the test without requiring assistance. 97.54% of the different types of results were interpreted correctly.

### LIMITATIONS

- The kit is a qualitative test that cannot quantify the concentration of SARS-CoV-2 antigen.
- The result of this test is not the only confirmation indicator for clinical indications. If the test result does not match clinical evidence, it is recommended to perform additional tests to verify the result.
- Children and young people under the age of 18 should be supported by an adult.
- This test cannot be used for non-human samples.

### ADDITIONAL REQUIRED EQUIPMENT

Need a clock or timer to see the time.

### PRECAUTIONS

- The kit is intended for in-vitro diagnostic use only. Please read the instructions for use carefully before test.
- Please use the swab and sample extraction buffer included in this kit. Do not replace the sample extraction in this kit with components from other kits.
- Use only undamaged test kits.
- Operation should be strictly according to the instructions and different batches should not be mixed.
- Do not drink the buffer at any time. If you drink, please consult a doctor immediately.

### EXPLANATION OF SYMBOL

Symbol	Meaning
Use by date	Consult instruction for use
Content Sufficient For re-tests	Temperature limitation
Manufacturing date	Caution
CE Marking - IVD 98/79/EC	Authorized representative in the European Community
For In Vitro Diagnostic Use	Keep away from sunlight
Do not use if package is damaged	Keep dry
CE mark	Sterilized using ethylene oxide

### ACCESSORIES

**What should I do if I did the test but don't see a control line?** In this case the test result is to be considered invalid. This can be caused by a possibly incorrect test operation. Please repeat the test with a new test kit according to the instructions for use. If you have any further invalid test results, contact your doctor or a COVID-19 test center.

**I am unsure of the interpretation of the results. What should I do?** If you cannot clearly determine the result of the test, contact the nearest medical facility using the rules of your local authority.

### My result is positive. What should I do?

If a horizontal colored line is visible in the control area (C) as well as in the test area (T), your result is positive. There is currently a suspicion of a COVID-19 infection. Immediately take self-isolation measures in accordance with the local guidelines and immediately contact your general practitioner / doctor or the local health department in accordance with the instructions of your local authorities. Your test result will be checked by a PCR confirmation test and the next steps will be explained to you.

### My result is negative. What should I do?

If only a horizontal colored line is visible in the control area (C), this can mean that you are negative or that the viral load is too low to be recognized by the test. If you experience symptoms such as headaches, migraines, fever, loss of sense of smell and taste, contact the nearest medical facility using the rules of your local authority. In addition, you can repeat the test with a new test kit. In case of suspicion, repeat the test after 1-2 days, as the coronavirus cannot be precisely detected in all phases of an infection. Distance and hygiene rules must still be observed.

### Can this test cassette be reused or used by more than one person?

This test cassette is for single use and cannot be reused or used by more than one person.

### Where do I dispose of the product?

The test kit can be disposed of with normal household waste.

Beijing Hotgen Biotech Co., Ltd.  
Wu Building, No. 9 Tianru Street, Biomedical Base, Daxing District, Beijing, 102602, P.R. China.

MedNet GmbH  
Birkfeldstraße 10, 48653 Münster, Germany

APPROVAL DATE AND REVISION DATE OF THE INSTRUCTION  
Approved on July 2021.  
Version number: 2021-07-01[Eng]

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